

SterilGARD



ENERGY EFFICIENT ENGINEERED

Class II, Type A2 Biological Safety Cabinet, Vertical Flow



SG403aHE
@ \$S10K
while-stock last

- Multiple energy saving features equals up to 60% increase in energy efficiency
- Larger and improved performance envelope for optimum safety
- Safest work environment with new StediFLOW™ self-adjusting motor technology
- With a 50% reduction in amperage draw and a quieter motor, noise is reduced up to 4 decibels vibration is minimized, and heat generation is reduced—making SterilGARD e3 the most comfortable working environment available
- An increase of up to 250% in the automatic motor reserve means 30% more filter and less impact on the environment
- More energy efficient and brighter lighting system provides 120 foot-candles for easier viewing and reduced eyestrain
- New, NSF-approved cable port in sidewall comes standard and allows for tubes and cables to be connected to equipment inside the work area without blocking the front grill
- The direct-to-plenum, filter-mounting technique with HEPA filters on a telescoping plenum makes maintenance easier and faster
- Industry's most durable and reliable cabinet means lower life cycle costs and years of trouble-free operation

For more information, please contact :

SciMed
the science of life.....

Tel : 6779 3388 Fax : 6266 3086

Email : sales@scimed.com.sg