

Systematic



- * Ergonomically Designed
- * Fully Autoclavable up to 121°C (250°F)
- * Color-Coded for Easy Identification on Type of Pipette Tips to Use
- * Designed to Fit Universal Pipette Tips
- * Easy Volume Adjustment for "A-Series" Models
- * Carousel Pipette Stand Available



Accessories

- * Pipette Rack - for 5 pipette (WH-5)
- * Pipette Rack - for 8 pipette (WH-8)



Precision Pipettes Selection Guide

A-Series - Variable Volume Precision Pipette

Model	Volume	Compatible Pipette Tip
A2	0.1-2.0uL	Clear
A10	1.0 to 10uL	Clear
A20	2.0 to 20uL	Yellow
A100	10 to 100uL	Yellow
A200	20 to 200uL	Yellow
A1000	100 to 1000uL	Blue
A5000	1000 to 5000uL	Clear

SGD 99

F-Series - Fixed Volume Precision Pipette

Model	Volume	Compatible Pipette Tip
F5	5uL	Clear
F10	10uL	Clear
F20	20uL	Yellow
F25	25uL	Yellow
F50	50uL	Yellow
F100	100uL	Yellow
F200	200uL	Yellow
F500	500uL	Blue
F1000	1000uL	Blue

SGD 95

* Each Pipette is complete with self-calibration kit

* Promotion valid whilst stock last



SciMed (Asia) Pte Ltd
196 Pandan Loop #07-11
Pantech Industrial Complex
Singapore 128384

Tel : 6779 3388 Fax: 6266 3086
Email : sales@scimed.com.sg

PRECISION MICRO PIPETTES